Contents :

α -Amylase	2
ALK-p	3
ALT	4
AST	5
B2M	6
Calcium	7
Cholestrol	8
CK-MB	9
Copper	10
CPK	11
Creatinine	12
CRP	13
D-Dimer	14
Direct Bilirubin	15
Glucose	16
HbA1C	17
HDL-C	19
Iron	20
LDH	21
LDL-C	22
Magnesium	23
Phosphorus	24
TIBC	25
Total Bilirubin	26
Triglycerides	27
Urea	28
Uric Acid	29
Urine Protein	30
Zinc	31

Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Amylase

Test	Amylase	R1	300
Reac.Type	Kinetic	R2	
Pri.Wave	412	Sample Volume	6
Sec.Wave		Linearity Range	0 2000
Direction	Increase	Factor	
Reac.Time	4 16	Slope	1
Incuba.Time		Intercept	0
Unit	U/L		
Precision	1		
Calibration	Rule	Two Point	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

ALP

Test	ALP	R1	200
Reac.Type	Kinetic	R2	50
Pri.Wave	450	Sample Volume	5
Sec.Wave		Linearity Range	0 1200
Direction	Increase	Factor	#
Reac.Time	4 16	Slope	1
Incuba.Time		Intercept	0
Unit	U/L		
Precision	1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

ALT

Test	ALT	R1	200
Reac.Type	Kinetic	R2	50
Pri.Wave	340	Sample Volume	25
Sec.Wave		Linearity Range	0 400
Direction	Decrease	Factor	
Reac.Time	4 16	Slope	1
Incuba.Time		Intercept	0
Unit	U/L		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

AST

Test	AST	R1	200
Reac.Type	Kinetic	R2	50
Pri.Wave	340	Sample Volume	25
Sec.Wave		Linearity Range	0 400
Direction	Decrease	Factor	
Reac.Time	4 16	Slope	1
Incuba.Time		Intercept	0
Unit	U/L		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Calcium

Test	Ca	R1	300
Reac.Type	End Point	R2	
Pri.Wave	670	Sample Volume	3
Sec.Wave		Linearity Range	0.5 16
Direction	Increase	Factor	
Reac.Time	0 18	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Cholesterol

Test	Chol	R1	300
Reac.Type	End point	R2	0
Pri.Wave	510	Sample Volume	3
Sec.Wave		Linearity Range	5 500
Direction	Increase	Factor	#
Reac.Time	0 36	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

CK-MB

Test	CK-MB	R1	200
Reac.Type	Kinetic	R2	50
Pri.Wave	340	Sample Volume	10
Sec.Wave		Linearity Range	0 318
Direction	Increase	Factor	#
Reac.Time	6 20	Slope	1
Incuba.Time	3	Intercept	0
Unit	U/L		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Copper

Test	Copper	R1	300
Reac.Type	End Point	R2	
Pri.Wave	578	Sample Volume	15
Sec.Wave	700	Linearity Range	12 400
Direction	Increase	Factor	
Reac.Time	0 18	Slope	1
Incuba.Time		Intercept	0
Unit	µg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	1	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

CPK

Test	CPK	R1	200
Reac.Type	Kinetic	R2	50
Pri.Wave	340	Sample Volume	7
Sec.Wave		Linearity Range	2 1000
Direction	Increase	Factor	#
Reac.Time	6 20	Slope	1
Incuba.Time	3	Intercept	0
Unit	U/L		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Creatinine

Test	Creat	R1	120
Reac.Type	Kinetic	R2	120
Pri.Wave	505	Sample Volume	25
Sec.Wave		Linearity Range	0.2 12.8
Direction	Increase	Factor	#
Reac.Time	7 14	Slope1	1
Inc.Time	4	Intercept	0
Unit	mg/dl		
Precision	0.1		
Calibration	Rule	Two Point	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

CRP

Test	CRP	R1	240
Reac.Type	Fixed Time	R2	60
Pri.Wave	578	Sample Volume	3
Sec.Wave		Linearity Range	2 80
Direction	Increase	Factor	#
Reac.Time	8 18	Slope	1
Inc.Time	9	Intercept	0
Unit	mg/L		
Precision	0.1		
Calibration	Rule	SPLINE (6 P)	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Direct Bilirubin

Test	D.B	R1	200
Reac.Type	End point	R2	50
Pri.Wave	546	Sample Volume	25
Sec.Wave		Linearity Range	0 15
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Incuba.Time	9	Intercept	0
Unit	mg /dl		
Precision	0.01		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Glucose

Test	Glu	R1	300
Reac.Type	End point	R2	0
Pri.Wave	510	Sample Volume	3
Sec.Wave		Linearity Range	5 400
Direction	Increase	Factor	#
Reac.Time	0 36	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Hb

Test	Hb	R1	180
Reac.Type	End point	R2	0
Pri.Wave	480	Sample Volume	12
Sec.Wave		Linearity Range	
Direction	Increase	Factor	#
Reac.Time	0 19	Slope	1
Incuba.Time		Intercept	0
Unit	mmol/l		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

A1C

$$\text{HbA1c\%} = (\text{A1c} / \text{Hb}) * 91.48 + 2.152$$

Test	A1C	R1	180
Reac.Type	End Point	R2	60
Pri.Wave	660	Sample Volume	12
Sec.Wave		Linearity Range	
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Inc.Time	18	Intersept	0
Unit	mmol/l		
Precision	0.01		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

HDL

Test	HDL	R1	240
Reac.Type	End Point	R2	80
Pri.Wave	578	Sample Volume	3
Sec.Wave		Linearity Range	15 150
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Inc.Time	9	Intersept	0
Unit	mg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Iron

Test	Iron	R1	200
Reac.Type	End Point	R2	50
Pri.Wave	578	Sample Volume	25
Sec.Wave	0	Linearity Range	5 400
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Incuba.Time	9	Intercept	0
Unit	µg /dl		
Precision	1		
Calibration	Rule	Two Point	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

LDH

Test	LDH	R1	200
Reac.Type	KINETIC	R2	50
Pri.Wave	340	Sample Volume	5
Sec.Wave		Linearity Range	10 1600
Direction	Decrease	Factor	#
Reac.Time	4 16	Slope	1
Incuba.Time	3	Intercept	0
Unit	U/L		
Precision	1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

LDL

Test	LDL	R1	240
Reac.Type	End Point	R2	80
Pri.Wave	578	Sample Volume	3
Sec.Wave		Linearity Range	0 500
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Incuba.Time	9	Intercept	0
Unit	mg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Magnesium

Test	Mg	R1	300
Reac.Type	End Point	R2	
Pri.Wave	546	Sample Volume	3
Sec.Wave		Linearity Range	0.5 5
Direction	Increase	Factor	
Reac.Time	0 18	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	0.01		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Phosphorus

Test	Phos	R1	300
Reac.Type	End Point	R2	
Pri.Wave	340	Sample Volume	3
Sec.Wave		Linearity Range	0.4 12.8
Direction	Increase	Factor	
Reac.Time	0 18	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

TIBC

Test	TIBC	R1	250
Reac.Type	End Point	R2	80
Pri.Wave	670	Sample Volume	18
Sec.Wave		Linearity Range	80 680
Direction	Decrease	Factor	#
Reac.Time	-1 18	Slope	1
Incuba.Time	9	Intercept	0
Unit	µg /dl		
Precision	1		
Calibration	Rule	Two Point	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Total Bilirubin

Test	T.B	R1	200
Reac.Type	End point	R2	50
Pri.Wave	546	Sample Volume	6
Sec.Wave		Linearity Range	0 20
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Incuba.Time	9	Intercept	0
Unit	mg /dl		
Precision	0.01		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Triglycerides

Test	TG	R1	300
Reac.Type	End point	R2	
Pri.Wave	510	Sample Volume	3
Sec.Wave		Linearity Range	5 700
Direction	Increase	Factor	#
Reac.Time	0 36	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Urea

Test	Urea	R1	240
Reac.Type	Kinetic	R2	60
Pri.Wave	340	Sample Volume	3
Sec.Wave		Linearity Range	5 160
Direction	Decrease	Factor	#
Reac.Time	7 14	Slope1	
Inc.Time	4	Intercept	0
Unit	mg/dl		
Precision	0.1		
Calibration	Rule	Two Point	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Uric Acid

Test	U.A	R1	240
Reac.Type	End Point	R2	60
Pri.Wave	546	Sample Volume	6
Sec.Wave		Linearity Range	1 14.4
Direction	Increase	Factor	#
Reac.Time	-1 18	Slope	1
Inc.Time	9	Intersept	0
Unit	mg /dl		
Precision	0.1		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

URINE AND CSF PROTEIN

Test	U.P	R1	300
Reac.Type	End Point	R2	
Pri.Wave	605	Sample Volume	10
Sec.Wave		Linearity Range	
Direction	Increase	Factor	
Reac.Time	0	Slope	1
Incuba.Time		Intercept	0
Unit	mg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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Persian Tajhiz System
Medical Equipment, Diagnostics and Consumables

MINDRAY BS 200-200E

Zinc

Test	Zinc	R1	300
Reac.Type	End Point	R2	
Pri.Wave	578	Sample Volume	19.5
Sec.Wave		Linearity Range	5 400
Direction	Increase	Factor	
Reac.Time	0 18	Slope	1
Incuba.Time		Intercept	0
Unit	µg /dl		
Precision	1.0		
Calibration	Rule	TWO POINT	
	Replicate	2	

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