

# **Zinc Calibrator**

Volume: 1 ml

Material: Liquid

## **Preface:**

Zinc Calibrator is a Liquid Calibrator which is prepared from human serum.

## **Preparation:**

Open the vials lid Carefully and be aware not to lose the contents of the vial.

The calibrator is ready to use.

Avoid shaking the vial violently and creating foam.

## **Storage and Stability:**

All the components of vial are stable until the expiration date on the label at 2-8°C.

## **Warning and precautions:**

This product is only for laboratory use.

**Add: NO.68, Sorena St., Damavand 2 Industrial Park, Damavand, Tehran, IRAN**

**Tell: 021-77778007-8**

**website: [www.pts-ico.ir](http://www.pts-ico.ir)**

**Fax: 021-26142195**

**email: [pts.ico@gmail.com](mailto:pts.ico@gmail.com)**

V.F1